

# Orbital Shaker

오비탈셰이커



## Features

- The height is low, little swing and no breakdown.
- Accurate speed control and digital mark function made for user convenience.
- Various holder and plate available according to users spec.
- Various flask usable by attaching Spring Wire Rack. (select)


## Specifications

Made in Korea

Spec. / Model		OS-300	OS-400
Motion		Orbital	
Speed	Range	25 to 300 rpm(up to 350rpm : high speed type)	
	Accuracy	±1%	
	Increment	1 rpm	
Time	Range	59min 59sec or Continuous	
	Accuracy	±1%	
	Increment	1 minute	
Controller		Digital PID Controller	
Display		4 Digit LED , 2 Display	
Drive system		Beltless Type	
Stroke		20mm	
Size	Plate (W x D mm)	330x330	470x470
	External (W x D x H, mm)	353x400x160	493x540x160
Electric Requirements		220VAC 50/60Hz	
Flask Number			
100mℓ		36ea	
250mℓ		23ea	
500mℓ		16ea	
1000mℓ		9ea	
2000mℓ		4ea	

## Obital Shaker

### Distinction



### Structure

