

# Water Bath

항온수조



## Features

- Equipped overheat protection sensor.
- Time set function and constant temperature.
- Uniquely designed smooth exterior and mono piece bath.
- Equipped internal bottom for the safety of the user and protection of the unit.

## Specifications


Made in Korea

Spec. / Model		WD-06	WD-11	WD-23
Type		Digital		
Controller		Digital PID Controller		
Range		Ambient +5~100 ° C		
Temp	Accuracy	±0.2 °C		
Size	Internal (W x D x H, mm)	240x140x150	300x240x150	515x330x150
	External (W x D x H, mm)	320x270x220	360x370x230	580x410x230
Safety		Over temp limiter		
Material	Internal	Seamless Stainless steel		
	External	Steel, Powder coating		
Electric Requirements		220VAC 50/60Hz		
Heater		500W	1000W	2400W

## Water Bath

항온수조

### Distinction



### Structure

