

Automatic Water Still

자동증류수기



Features

- Automatic Water Still produces and reserves distilled water automatically.
- With its fine design and jumbo size reservoir, this unit offers exceptional value for the user.
- It is very easy and safe to operate with indicators such as reservoir full indicator, distilling indicator, and with safety device such as overheat safety sensor.

Specifications


Made in Korea

Model		HWS-FA12	HWS-FA30
Capacity		3.5ℓ~4.0ℓ/ hour	3.5ℓ~4.0ℓ/ hour
Dimensions	Boiling Water Tank (W x D x H, mm)	200 x 200 x 300	200 x 200 x 300
	Distilled Water Tank (W x D x H, mm)	200 x 200 x 300	275 x 275 x 400
	Out (W x D x H, mm)	530 x 300 x 590	680 x 375 x 600
Heater		3.0 kW	
Distilled Water		Approx. 8ℓ	Approx. 20ℓ
Material	Internal	Stainless Steel (#304)	
	External	Powder Coated Steel	
Safety Device		Over Temp Limiter	
		Water Supply Cut-off	
Operation		Fully Automatic Operation	
Signal Lamp		Distillation, Filled with water Lamp	
Electric Supply		220V, 50/60 Hz	

Automatic Water Still

자동증류수기

Distinction



Structure

