

구매 규격서

Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		High Temp Oven 고온건조기 (C-DHD2)	Set	1

1. 용도 및 특징 (Features)

시료,초자기구등의 세척 후 건조등 다방면으로 사용되는 건조기입니다.

- 자연대류형 건조기로 고온의 건조기능이 필요할 때 사용되는 제품입니다.
- 작은크기의 실험물, 가루 등 바람에 영향을 받지 않고 고온에 시료를 건조 시키기 유리합니다.
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있습니다.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗합니다.
- 온도 Safety 장치, 과전류, 누전 장치등 안전장치 강화

2. 사양 (Specifications)

Model		C-DHD2
Dimensions	Inner (W*D*H)(mm)	400×400×450
	Outer (W*D*H)(mm)	520×540×850
	Volume (L)	72
Controller		Digital Temp Controller & Digital Display / K type(CA)
Temp	Range	Ambient +5℃ ~ 300 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	2,000W
	Power	220V, 9A, 1ø
Shelf		3EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		58

3. 구성 (Composition)

- Manual..... 1ea
- Shelf..... 1set
- Power cord cable..... 1ea

4. 선택옵션 (Optional Accessories)

- Perforated shelves
- Basket shelves

5. 보증 (Warranty)

- 1 Year

구매 규격서

Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		High Temp Oven 고온건조기 (C-DHD3)	Set	1

1. 용도 및 특징 (Features)

시료,초자기구등의 세척 후 건조등 다방면으로 사용되는 건조기입니다.

- 자연대류형 건조기로 고온의 건조기능이 필요할 때 사용되는 제품입니다.
- 작은크기의 실험물, 가루 등 바람에 영향을 받지 않고 고온에 시료를 건조 시키기 유리합니다.
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있습니다.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗합니다.
- 온도 Safety 장치, 과전류, 누전 장치등 안전장치 강화

2. 사양 (Specifications)

Model		C-DHD3
Dimensions	Inner (W*D*H)(mm)	500×450×600
	Outer (W*D*H)(mm)	620×590×1000
	Volume (L)	135
Controller		Digital Temp Controller & Digital Display / K type(CA)
Temp	Range	Ambient +5℃ ~ 300 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	3,000W
	Power	220V, 13.6A, 1ø
Shelf		3EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		80

3. 구성 (Composition)

- Manual..... 1ea
- Shelf..... 1set
- Power cord cable..... 1ea

4. 선택옵션 (Optional Accessories)

- Perforated shelves
- Basket shelves

5. 보증 (Warranty)

- 1 Year