

구매 규격서

Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Shaking Incubator 항온진탕배양기 (C-SKI-1)	Set	1

1. 용도 및 특징 (Features)

- 회전형으로(ORBITAL방식) Shaking , 0℃ ~60℃ 까지 사용
- Digital Temp Controller & Digital R.P.M Controller 로 정확한 온도, 속도 조절 가능
- AC Motor 를 사용 저소음으로 장기간 사용시 에도 내구성이 뛰어납니다.
- Feed-Back Speed Controller를 채용하여 전압, 부하의 변동에 변화되지 않은 R.P.M을 유지
- 과부하차단장치, 누전차단장치, Door Open시 자동정지 등 안전장치가 내장
- 500ml Flask Holder × 12ea 기본 장착 (Holder size 변경 가능, Spring rack 교체가능)

2. 사양 (Specifications)

Model		C-SKI-1
Dimensions	Inner (W*D*H)(mm)	540×540×400
	Outer (W*D*H)(mm)	590×780×930
Controller		Digital Temp Controller & Digital Speed controller / PT 100Ω
R.P.M		0 ~ 300
Temp	Range	0 ℃ ~ 60 ℃
	Accuracy	Better then ±0.2 at 37℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Flask Holder		500ml × 12ea
Electrical spec	Power	220V, 300W , 1ø
Materials	Plate	Stainless steel
	Body	Steel plate with powder coated
	Door	Acrylic
Refrigerator		1/4 HP
Safety device		Exclusive over temp, protector , Door opening auto stop
Weignt(Kg)		90

3. 구성 (Composition)

- Manual..... 1ea
- Flask Holder..... 1set

4. 선택옵션 (Optional Accessories)

- Flask Holder
- Spring Rack
- A lighting installation

5. 보증 (Warranty)

- 1 Year

구매 규격서

Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Shaking Incubator 항온진탕배양기 (C-SKI-2)	Set	1

1. 용도 및 특징 (Features)

- 회전형으로(ORBITAL방식) Shaking , Ambient+5℃ ~ 60℃ 까지 사용
- Digital Temp Controller & Digital R.P.M Controller 로 정확한 온도, 속도 조절 가능
- AC Motor 를 사용 저소음으로 장기간 사용시 예도 내구성이 뛰어납니다.
- Feed-Back Speed Controller를 채용하여 전압, 부하의 변동에 변화되지 않은 R.P.M을 유지
- 과부하차단장치, 누전차단장치, Door Open시 자동정지 등 안전장치가 내장
- 500ml Flask Holder × 12ea 기본 장착 (Holder size 변경 가능, Spring rack 교체가능)

2. 사양 (Specifications)

Model		C-SKI-2
Dimensions	Inner (W*D*H)(mm)	500×480×260
	Outer (W*D*H)(mm)	530×800×520
Controller		Digital Temp Controller & Digital Speed controller / PT 100Ω
R.P.M		0 ~ 300
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	Better then ±0.2 at 37℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Flask Holder		500ml × 12ea
Electrical spec	Power	220V, 300W , 1ø
Materials	Plate	Stainless steel
	Body	Steel plate with powder coated
	Door	Acrylic
Safety device		Exclusive over temp, protector , Door opening auto stop
Weignt(Kg)		50

3. 구성 (Composition)

- Manual..... 1ea
- Flask Holder..... 1set

4. 선택옵션 (Optional Accessories)

- Flask Holder
- Spring Rack
- A lighting installation

5. 보증 (Warranty)

- 1 Year