

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Vacuum Oven 진공건조기 (C-DVD1)	Set	1

### 1. 용도 및 특징 (Features)

식품의 건조 및 수분 측정시 수분제거용, 분말제품의 탈기, 용매제거, 살균건조용입니다.

- 직접가열이 어려운 시료의 건조에 사용됩니다.
- 내부재질은 SUS304 Plate 3t를 사용하고 내부는 깨끗이 마감
- 외부재질은 철 판에 분체코팅 하여 제품이 견고하고 미려함
- Aluminum shelf 은 열전도가 우수
- 내통의 외부를 감싸 간접으로 가열하는 방식 제작되어 내부의 온도분포가 균일
- 겉문의 특수제작된 실리콘패킹으로 진공도가 우수
- 진공게이지 부착으로 실험 시 진공도를 확인 가능
- 누전차단기를 부착하여 과전류와 누전이 될 때 전원이 차단

### 2. 사양 (Specifications)

Model		C-DVD1
Dimensions	Inner (W*D*H)(mm)	300×300×350
	Outer (W*D*H)(mm)	450×410×630
	Volume (L)	27
Controller		Digital Temp Controller & Digital Display / K type(CA)
Temp	Range	Ambient +5℃ ~ 250 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	800W
	Power	220V, 3.6A, 1ø
Shelf		1EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		48

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelf..... 1ea
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Aluminium shelf
- Vacuum Pump
- Cold Trap bath

### **5. 보증 (Warranty)**

- 1 Year

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Vacuum Oven 진공건조기 (C-DVD2)	Set	1

### 1. 용도 및 특징 (Features)

식품의 건조 및 수분 측정시 수분제거용, 분말제품의 탈기, 용매제거, 살균건조용입니다.

- 직접가열이 어려운 시료의 건조에 사용됩니다.
- 내부재질은 SUS304 Plate 3t를 사용하고 내부는 깨끗이 마감
- 외부재질은 철판에 분체코팅 하여 제품이 견고하고 미려함
- Aluminum shelf 은 열전도가 우수
- 내통의 외부를 감싸 간접으로 가열하는 방식 제작되어 내부의 온도분포가 균일
- 겉문의 특수제작된 실리콘패킹으로 진공도가 우수
- 진공게이지 부착으로 실험 시 진공도를 확인 가능
- 누전차단기를 부착하여 과전류와 누전이 될 때 전원이 차단

### 2. 사양 (Specifications)

Model		C-DVD2
Dimensions	Inner (W*D*H)(mm)	350×350×350
	Outer (W*D*H)(mm)	830×525×670
	Volume (L)	43
Controller		Digital Temp Controller & Digital Display / K type(CA)
Temp	Range	Ambient +5℃ ~ 250 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	1000W
	Power	220V, 4.5A, 1ø
Shelf		Aluminium Plate Shelf / 1EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		58

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelf..... 1ea
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Aluminium shelf
- Vacuum Pump
- Cold Trap bath

### **5. 보증 (Warranty)**

- 1 Year

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Vacuum Oven 진공건조기 (C-DVD3)	Set	1

### 1. 용도 및 특징 (Features)

식품의 건조 및 수분 측정시 수분제거용, 분말제품의 탈기, 용매제거, 살균건조용입니다.

- 직접가열이 어려운 시료의 건조에 사용됩니다.
- 내부재질은 SUS304 Plate 3t를 사용하고 내부는 깨끗이 마감
- 외부재질은 철판에 분체코팅 하여 제품이 견고하고 미려함
- Aluminum shelf 은 열전도가 우수
- 내통의 외부를 감싸 간접으로 가열하는 방식 제작되어 내부의 온도분포가 균일
- 겉문의 특수제작된 실리콘패킹으로 진공도가 우수
- 진공게이지 부착으로 실험 시 진공도를 확인 가능
- 누전차단기를 부착하여 과전류와 누전이 될 때 전원이 차단

### 2. 사양 (Specifications)

Model		C-DVD3
Dimensions	Inner (W*D*H)(mm)	400×400×400
	Outer (W*D*H)(mm)	880×575×720
	Volume (L)	64
Controller		Digital Temp Controller & Digital Display / K type(CA)
Temp	Range	Ambient +5℃ ~ 250 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	1,200W
	Power	220V, 5.4A, 1ø
Shelf		Aluminium Plate Shelf / 1EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		74

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelf..... 1ea
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Aluminium shelf
- Vacuum Pump
- Cold Trap bath

### **5. 보증 (Warranty)**

- 1 Year