

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Normal Incubator 일반항온기 (C-IND1)	Set	1

### 1. 용도 및 특징 (Features)

열과 광선에 민감한 배양생물의 보존 연구실험용, 항생물질, 효소등의 배양, 균배양, 발아 시험용으로 쓰임

- 자연대류형 항온기로 일반적인 온도외에 발아 균을 배양
- 이중 문으로 제작, 겉문을 열고 속문을 열지 않은 상태에서 내부를 관찰 할 때 내부의 온도가 떨어지는것을 방지
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있음.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작이 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗 합니다.

### 2. 사양 (Specifications)

Model		C-IND1
Dimensions	Inner (W*D*H)(mm)	300×300×350
	Outer (W*D*H)(mm)	430×400×600
	Volume (L)	30
Controller		Digital Temp Controller & Digital Display / PT 100Ω
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	150W
	Power	220V, 0.7A, 1ø
Shelves		2EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		23.5

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelves..... 1set
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Perforated shelves
- Basket shelves

### **5. 보증 (Warranty)**

- 1 Year

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Normal Incubator 일반항온기 (C-IND2)	Set	1

### 1. 용도 및 특징 (Features)

열과 광선에 민감한 배양생물의 보존 연구실험용, 항생물질, 효소등의 배양, 균배양, 발아 시험용으로 쓰임

- 자연대류형 항온기로 일반적인 온도외에 발아 균을 배양
- 이중 문으로 제작, 겉문을 열고 속문을 열지 않은 상태에서 내부를 관찰 할 때 내부의 온도가 떨어지는것을 방지
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있음.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작이 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗 합니다.

### 2. 사양 (Specifications)

Model		C-IND2
Dimensions	Inner (W*D*H)(mm)	400×400×450
	Outer (W*D*H)(mm)	500×500×600
	Volume (L)	72
Controller		Digital Temp Controller & Digital Display / PT 100Ω
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	250W
	Power	220V, 1A, 1ø
Shelves		3EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		35

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelves..... 1set
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Perforated shelves
- Basket shelves

### **5. 보증 (Warranty)**

- 1 Year

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Normal Incubator 일반항온기 (C-IND3)	Set	1

### 1. 용도 및 특징 (Features)

열과 광선에 민감한 배양생물의 보존 연구실험용, 항생물질, 효소등의 배양, 균배양, 발아 시험용으로 쓰임

- 자연대류형 항온기로 일반적인 온도외에 발아 균을 배양
- 이중 문으로 제작, 겉문을 열고 속문을 열지 않은 상태에서 내부를 관찰 할 때 내부의 온도가 떨어지는것을 방지
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있음.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작이 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗 합니다.

### 2. 사양 (Specifications)

Model		C-IND3
Dimensions	Inner (W*D*H)(mm)	500×500×600
	Outer (W*D*H)(mm)	590×550×850
	Volume (L)	150
Controller		Digital Temp Controller & Digital Display / PT 100Ω
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	±1℃
Timer	Scale	99Hr 59Min, 99Min 59Sec, 99Day 23Hr
Electrical spec	Heater	250W
	Power	220V, 1A, 1ø
Shelves		3EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		40

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelves..... 1set
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Perforated shelves
- Basket shelves

### **5. 보증 (Warranty)**

- 1 Year

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Normal Incubator 일반항온기 (C-INA1)	Set	1

### 1. 용도 및 특징 (Features)

열과 광선에 민감한 배양생물의 보존 연구실험용, 항생물질, 효소등의 배양, 균배양, 발아 시험용으로 쓰임

- 자연대류형 항온기로 일반적인 온도외에 발아 균을 배양
- 이중 문으로 제작, 걸문을 열고 속문을 열지 않은 상태에서 내부를 관찰 할 때 내부의 온도가 떨어지는것을 방지
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있음.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작이 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗 합니다.

### 2. 사양 (Specifications)

Model		C-INA1
Dimensions	Inner (W*D*H)(mm)	300×300×350
	Outer (W*D*H)(mm)	430×400×600
	Volume (L)	30
Controller		Analog Temp Controller
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	±5℃
Timer	Scale	×
Electrical spec	Heater	150W
	Power	220V, 0.7A, 1ø
Shelves		2EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		23.5

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelves..... 1set
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Perforated shelves
- Basket shelves

### **5. 보증 (Warranty)**

- 1 Year



# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Normal Incubator 일반항온기 (C-INA2)	Set	1

### 1. 용도 및 특징 (Features)

열과 광선에 민감한 배양생물의 보존 연구실험용, 항생물질, 효소등의 배양, 균배양, 발아 시험용으로 쓰임

- 자연대류형 항온기로 일반적인 온도외에 발아 균을 배양
- 이중 문으로 제작, 걸문을 열고 속문을 열지 않은 상태에서 내부를 관찰 할 때 내부의 온도가 떨어지는것을 방지
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있음.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작이 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗 합니다.

### 2. 사양 (Specifications)

Model		C-INA2
Dimensions	Inner (W*D*H)(mm)	400×400×450
	Outer (W*D*H)(mm)	500×500×600
	Volume (L)	72
Controller		Analog Temp Controller
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	±5℃
Timer	Scale	×
Electrical spec	Heater	250W
	Power	220V, 1A, 1ø
Shelves		3EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		35

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelves..... 1set
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Perforated shelves
- Basket shelves

### **5. 보증 (Warranty)**

- 1 Year

# 구매 규격서

## Commodity Description

관세 분류번호 H. S. K NO	품목 번호 Item. NO	품명 및 규격 Description	단위 Unit	수량 Q'TY
		Normal Incubator 일반항온기 (C-INA3)	Set	1

### 1. 용도 및 특징 (Features)

열과 광선에 민감한 배양생물의 보존 연구실험용, 항생물질, 효소등의 배양, 균배양, 발아 시험용으로 쓰임

- 자연대류형 항온기로 일반적인 온도외에 발아 균을 배양
- 이중 문으로 제작, 겔문을 열고 속문을 열지 않은 상태에서 내부를 관찰 할 때 내부의 온도가 떨어지는것을 방지
- 내부를 SUS광판(光板)을 사용하여 깨끗하게 마감이 되어 있음.
- 외부는 철판에 분체코팅을 하여 제품이 견고하고 미려하며 방청효과가 우수
- 선반이 WIRE로 제작이 되어 내부의 공기가 순환하기 원활하게 제작
- 실리콘 패킹을 사용 밀폐성이 좋고 깨끗

### 2. 사양 (Specifications)

Model		C-INA3
Dimensions	Inner (W*D*H)(mm)	500×500×600
	Outer (W*D*H)(mm)	590×550×850
	Volume (L)	150
Controller		Analog Temp Controller
Temp	Range	Ambient +5℃ ~ 60 ℃
	Accuracy	±5℃
Timer	Scale	×
Electrical spec	Heater	250W
	Power	220V, 1A, 1ø
Shelves		3EA
Materials	Inner	Stainless steel SUS304
	Outer	Steel plate with powder coated
Safety device		Exclusive over temp, protector
Weignt(Kg)		40

### **3. 구성 (Composition)**

- Manual..... 1ea
- Shelves..... 1set
- Power cord cable..... 1ea

### **4. 선택옵션 (Optional Accessories)**

- Perforated shelves
- Basket shelves

### **5. 보증 (Warranty)**

- 1 Year